

NEW

ASPIRATION RANGE


CORAIL PLUS

BALLOON-TIPPED GUIDING CATHETER 6F AND 8F



VASCO PLUS

BRAIDED MICROCATETER



FARGO



L'INNOVATION SUR MESURE*

www.balt.fr


ASPIRATION RANGE

BALT is pleased to present its new range of thrombus aspiration devices specifically designed for different clot locations.

CORAIL8F+ASPI

Balloon-tipped 8F guiding catheter:

- ✓ For proximal aspiration in the **Carotid Artery**
- ✓ Distal balloon for **temporary occlusion**
- ✓ Large inner lumen allowing **easy navigation for a catheter and compatible with all thrombectomy devices**
- ✓ Reference: CORAIL8F+




VASCO+35ASPI

Braided microcatheter :

- ✓ For distal aspiration up to the **Middle Cerebral Artery**
- ✓ **Progressive suppleness and extra supple tip**
- ✓ **Large inner lumen (1.02mm, 3.1F)**
- ✓ Profile especially adjusted to **limit the aspiration of the upstream blood flow**
- ✓ Reference: VASCO+35ASPI (D/MP)

I.D.:1.02mm (.0401")




VASCO+35ASPI can be inserted into the CORAIL8F+ASPI for a combined aspiration effect

FARGO6F ASPI

6F guiding catheter:

- ✓ For thrombus aspiration in the **Posterior Circulation**
- ✓ **Distal outer diameter (4.9F)** for flow control in the basilar artery
- ✓ **Large inner lumen (1.40mm, 4.2F)**
- ✓ **Long ultra-supple and navigable 15 CM** distal tip
- ✓ Reference: FRG6F135_15 (D/MP)



-Use a long introducer IVA6F80 to position the FARGO6F ASPI
-FARGO6F ASPI can also be used with a VASCO+28 or any other catheter with an O.D of 3.4F



10, rue Croix Vigneron
95160 Montmorency - France
Tél : +33 1 39 89 46 41
Fax : +33 1 34 17 03 46
www.balt.fr