

DIALYSIS

DIALYSIS SET

Single lumen, straight	2
Single lumen "Y"	2
Double lumen	3
Triple lumen	4

DIALYSIS CATHETER

Single lumen	5
Single lumen "Y"	5
Double lumen	6
Triple lumen	6

HEMODIALYSIS NEEDLE & ACCESSORIES

Hemodialysis needle	7
Dialysis extender for dialyzer rinsing / cleaning	8
Dialysis extender "Y"	8

DIALYSIS SET single lumen, straight



Set components:

single lumen straight catheter
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size / length	Dilator I / II	Guide wire size / length
ZDJP8F15	ZDJP8F15P	8F / 15 cm	8F /10F	.035" / 60 cm
ZDJP8F18	ZDJP8F18P	8F / 18 cm	8F /10F	.035" / 60 cm
ZDJP8F20	ZDJP8F20P	8F / 20 cm	8F /10F	.035" / 60 cm

When ordering use Ref. Number, e.g.: ZDJP8F20.

NOTE: On request we provide catheters with hydrophilic surface protection. When ordering use Ref. Number, e.g.: ZDJP8F20H.

DIALYSIS SET single lumen "Y"



Set components:

single lumen catheter "Y"
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size / length	Dilator I / II	Guide wire size / length
ZDZY8F15	ZDZY8F15P	8F / 15 cm	8F /10F	.035" / 60 cm
ZDZY8F18	ZDZY8F18P	8F / 18 cm	8F /10F	.035" / 60 cm
ZDZY8F20	ZDZY8F20P	8F / 20 cm	8F /10F	.035" / 60 cm


When ordering use Ref. Number, e.g.: ZDZY8F20.

NOTE: On request we provide catheters with hydrophilic surface protection. When ordering use Ref. Number, e.g.: ZDZY8F20H.

DIALYSIS SET double lumen

Set components:

double lumen catheter Straight type
 double lumen catheter Curved type
 needle 18G, length 7 cm
 for catheter size 6.5F - needle 19G, length 5 cm
 guide wire
 dilatator I
 dilatator II



Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size / length	Dilatator I / II	Guide wire size / length	Version
ZDD6.5F10	ZDD6.5F10P	6.5F / 10 cm	7F / 8F	.027" / 60 cm	PEDIATRIC
ZDD6.5F12	ZDD6.5F12P	6.5F / 12 cm	7F / 8F	.027" / 60 cm	
ZDD8F10	ZDD8F10P	8F / 10 cm	8F / 10F	.035" / 60 cm	
ZDD8F12	ZDD8F12P	8F / 12 cm	8F / 10F	.035" / 60 cm	
ZDD8F15	ZDD8F15P	8F / 15 cm	8F / 10F	.035" / 60 cm	
ZDD11F15	ZDD11F15P	11F / 15 cm	10F / 12F	.035" / 70 cm	
ZDD11F18	ZDD11F18P	11F / 18 cm	10F / 12F	.035" / 70 cm	
ZDD11F20	ZDD11F20P	11F / 20 cm	10F / 12F	.035" / 70 cm	
ZDD12F15	ZDD12F15P	12F / 15 cm	10F / 12F	.035" / 70 cm	
ZDD12F18	ZDD12F18P	12F / 18 cm	10F / 12F	.035" / 70 cm	
ZDD12F20	ZDD12F20P	12F / 20 cm	10F / 12F	.035" / 70 cm	

When ordering use Ref. Number, e.g.: ZDD12F20.

NOTE: other configurations available


H - catheter with hydrophilic coating

When ordering use Ref. Number, e.g.: ZDD11F20HP.

Z - catheter with curved arms

When ordering use Ref. Number, e.g.: ZDD11F20ZP.

DIALYSIS SET triple lumen



Set components:

triple lumen catheter Straight type
 triple lumen catheter Curved type
 needle 18G, length 7 cm
 guide wire
 dilator I
 dilator II

Reference:	Catheter size / length	Dilatator I / II	Guide wire size / length
ZDT12F15	12F / 15 cm	10F / 12F	.038" / 70 cm
ZDT12F18	12F / 18 cm	10F / 12F	.038" / 70 cm
ZDT12F20	12F / 20 cm	10F / 12F	.038" / 70 cm

When ordering use Ref. Number, e.g.: ZDT12F15.

NOTE: other configurations available:

H - catheter with hydrophilic coating

When ordering use Ref. Number, e.g.: ZDT12F15H

Z - catheter with curved arms

When ordering use Ref. Number, e.g.: ZDT12F15Z.

DIALYSIS CATHETER single lumen



Reference:

KDJP8F15
KDJP8F18
KDJP8F20

Catheter size / length

8F / 15 cm
8F / 18 cm
8F / 20 cm

When ordering use Ref. Number, e.g.: KDJP8F20.

NOTE: Catheter with curved arms available to order. When ordering use Ref. Number, e.g.: KDJP8F15Z.

DIALYSIS CATHETER single lumen "Y"



Reference:


KDJY8F15
KDJY8F18
KDJY8F20

Catheter size / length

8F / 15 cm
8F / 18 cm
8F / 20 cm

When ordering use Ref. Number, e.g.: KDJY8F20.

DIALYSIS CATHETER double lumen



Reference:	Catheter size / length	Version
KDD6.5F10	6.5F / 10 cm	PEDIATRIC
KDD6.5F12	6.5F / 12 cm	
KDD8F10	8F / 10 cm	
KDD8F12	8F / 12 cm	
KDD8F15	8F / 15 cm	
KDD11F15	11F / 15 cm	FOR ADULTS
KDD11F18	11F / 18 cm	
KDD11F20	11F / 20 cm	
KDD12F15	12F / 15 cm	
KDD12F18	12F / 18 cm	
KDD12F20	12F / 20 cm	

When ordering use Ref. Number, e.g.: KDD8F12.

NOTE: Catheter with curved arms available to order. When ordering use Ref. Number, e.g.: KDD8F15Z.
H - catheter with hydrophilic coating. When ordering use Ref. Number, e.g.: KDD12F15H.

DIALYSIS CATHETER triple lumen



Reference:	Catheter size / length
KDT12F15	12F / 15 cm
KDT12F18	12F / 18 cm
KDT12F20	12F / 20 cm

When ordering use Ref. Number, e.g.: KDT12F15.

NOTE: other configurations available
H - catheter with hydrophilic coating. When ordering use Ref. Number, e.g.: KDT12F15H
Z - catheter with curved arms. When ordering use Ref. Number, e.g.: KDT12F15Z.

HEMODIALYSIS NEEDLE

The needle is made of special steel alloy that ensures exceptional smoothness of its surface. Thin walls and a triple - grinded spike punctures the fistula delicately and precisely. The surface of the needle is covered with silicone and therefore it goes through the tissues easily and well tolerated in the vessel. Specially shaped wings hold the needle easily and allow easy hanging. Moving wings make it easy to attach the needle to the skin - the position of the butterfly - nut can be changed without turning the needle in the vessel. The clamps that are easy to fasten and reliable in action are colored according to the line concerned: blue in the case of venous needles, red in the case of arterial needles, that allows easy application and prevents mistakes.



Reference:	Needle size	Version	Wings colour	Clamps colour
ITD1.5LL	17G	Arterial	red	red
ITD1.6LL	16G	Arterial	green	red
ITD1.8LL	15G	Arterial	yellow	red
IZD1.5LL	17G	Venous	red	blue
IZD1.6LL	16G	Venous	green	blue
IZD1.8LL	15G	Venous	yellow	blue
IDY1.5LL	17G	"Y"	red	red/blue
IDY1.6LL	16G	"Y"	green	red/blue
IDY1.8LL	15G	"Y"	yellow	red/blue

When ordering use Ref. Number, e.g.: IZD1.5LL.

DIALYSIS EXTENDER for dializator rinsing / cleaning


Reference:

PD7

Length


7 cm

Diameter

4.6 / 3.3 mm

When ordering use Ref. Number, e.g.: PD7.

DIALYSIS EXTENDER "Y"



"Y" extender is used to transform straight catheter into a "Y" catheter.

Reference:

PUDY

When ordering use Ref. Number, e.g.: PUDY.